Manual de uso	ES
User manual	EN
Manuale d'uso	IT
Manual do utilizador	PT
Manual do utilizador	BR
Benutzerhandbuch	DE
Manuel de l'utilisateur	FR
Podręcznik użytkownika	PL
Felhasználói kézikönyv	HU
Manual de utilizare	RO
Användarmanual	SV
Brukerhåndbok	NO
Naudotojo vadovas	LT
Ръководство за потребителя	BG
Používateľská príručka	SK
Käyttöohje	FI
Упутство за употребу	SR
Kasutusjuhend	ET
Kullanım Kılavuzu	TR

IRUDEK 2000 S.L.
Pol. Erribera 8A
20150 Aduna (Guipúzcoa)
España
Tfno: +34 943 69 26 17
Fax: +34 943 69 25 26
irudek@irudek.com

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EN

Read the operating instructions carefully before using the stretcher and the rescue loop, train yourself properly, familiarise yourself with it and use it responsibly. Activities at height involve serious risks not outlined in this manual, where each user is responsible for the management of such risks, their safety, their actions and the consequences of these, if you do not assume this or do not understand this manual, do not use the equipment.

MEDICAL CONDITIONS THAT MAY AFFECT THE SAFETY OF THE USER:
UNDER NORMAL CONDITIONS OF USE, CARDIOVASCULAR PROBLEMS, RESPIRATORY DISEASES, MUSCULOSKELETAL DISORDERS AFFECTING THE SPINE, HIPS OR KNEES, OBESITY OR EXCESS WEIGHT, AS WELL AS NEUROLOGICAL OR BALANCE DISORDERS AND PERIPHERAL CIRCULATORY PROBLEMS HINDERING VEIN RETURN MUST BE TAKEN INTO ACCOUNT. IN AN EMERGENCY SITUATION FOLLOWING THE ARREST OF A FALL, THE RISKS ARE PARTICULARLY SERIOUS IN PEOPLE PREDISPOSED TO HARNESSE SYNDROME, WITH A HISTORY OF TRAUMA OR INJURY, WITH COAGULATION PROBLEMS OR ON ANTICOAGULANT TREATMENT, AND IN THOSE WITH NEUROLOGICAL CONDITIONS THAT CAN CAUSE SEIZURES.

DESCRIPTION

Flexible stretcher made of PVC combined with polyester.

The EMAN ROLL is a roll-up stretcher that allows the casualty to be transported both horizontally and vertically.

For vertical manoeuvres, the casualty must have a fall arrest system independent of the stretcher.

The stretcher incorporates an EN 1498 Type A rescue lanyard, which must be used in conjunction with a fall arrest system, except in situations where the use of a separate fall arrest system would present a greater risk than the use of the EN 1498 Type A rescue lanyard alone.

The stretchers are made of highly durable and flexible plastic and have dimensions of 200 x 90 cm. This plastic material is distinguished by its high abrasion resistance, does not discolor with physiological fluids, has a surface with antibacterial properties and withstands exposure to UV radiation.

All equipment is stored in an ergonomic backpack specially designed for easy transport, both outdoors and in confined spaces.

Content of the EMAN ROLL set:

4 additional handles to allow more people to assist in lifting the casualty, 4 straps for closing the stretcher, 2 straps for restraint, casualty, stretcher, straps for closing the head and foot straps, perimeter rope, pull strap, 2 suspension straps, rucksack.

Information on this equipment is available at the following link: <http://www.irudek.com>

NOMENCLATURE

1. Handles, 2. Straps for closing the stretcher, 3. Straps for securing the injured person, 4. Stretcher, 5. Straps for closing the head and foot sections, 6. Perimeter rope, 7. Drag straps, 8. Suspension straps, 9. Backpack, (10, 11 and 12 see Rescue Loops EN 1498)

ASSEMBLY

Fig.1: Remove the stretcher from the backpack.

Fig.2: Unroll the stretcher and roll both ends of the stretcher to the side opposite to where it was stored, until the stretcher can lie flat on the floor.

Fig.3: If necessary finish shaping the spikes on the ground, if you are going to slide into an obstacle, the spikes cannot be an obstacle.

Fig.4 a and b: If the locking straps are not in position, thread them through the flat holes and attach them to the stretcher with a lar's foot knot.

Fig.5: Disconnect the buckles and bring the straps to the sides of the stretcher to leave the stretcher space free.

Fig.6: When placing the casualty on the stretcher, choose the most appropriate technique for the situation. In this manual we will not describe the different techniques that exist for the correct handling and transfer of casualties.

Fig.7: Connect the 2 straps for securing the casualty.

Fig.8: Connect the 4 straps for the head and foot end closure by lifting the sides of the stretcher.

Fig.9: Connect the foot straps.

Fig.10: For horizontal displacements by dragging the stretcher, traction of the traction belts.

Fig.11: For horizontal movements, when lifting the stretcher by hand, use the four additional handles.

Fig.12: For suspended and vertical movements, use the anchor point of the rope at the head of the stretcher. To position the rope, find the midpoint and tie an 8 knot. The victim must wear a fall arrest system independent of the stretcher's progression.

Fig.13: If the stretcher is to be used in an upright position, it is advisable to fold the upper part of the stretcher with the help of the traction straps, in order to avoid knocks and impacts caused by falling objects.

Fig.14: Method 1: Pass the ends through the eyelets at the head of the stretcher from the outside to the inside of the stretcher.

Fig.15: Method 2 for extremely small openings, start by passing the end of the rope inside out from the shoulder eyelets. The knot position should be the same as in Method 1.

Fig.16: Pass the rope through all the stretcher eyelets up to the foot slots.

Fig.17: Join the two ends of the rope with a flat knot in the middle of the stretcher.

Fig.18: Pass the rope through the lower handles and tie a re-secured flat knot in the middle of the stretcher.

Fig.19: Example of EMAN ROLL stretcher suspension progression in vertical position, including fall arrest system independent of the progression.

Fig.20 a and b: For suspended and horizontal movements, place the suspension straps through the eyelets at the waist and weight of the stretcher, as shown in the picture. Long strap at the feet, short strap at the waist. Avoid twisting, keeping the strap flat all the way around.

Fig.21: Use the anchorage point formed by the lugs of the suspension straps. It is recommended that a multi-anchorage plate be used to connect the two karabiners of the lanyards. The casualty must wear a fall arrest system independent of the stretcher's progression.

Fig.22 a and b: For suspended transfers from horizontal to vertical and vice versa, you can use the EMAN Stif or EMAN Stif Pro device, which is not included with the stretcher.

LIMITATIONS ON USE

The device may only be handled by appropriately trained personnel.

THE STRETCHER IS NOT A FALL ARREST SYSTEM AND THEREFORE, WHEN PERFORMING MANOEUVRES WITH A RISK OF FALLING TO A DIFFERENT LEVEL, THE PERSON OCCUPYING THE STRETCHER MUST BE CORRECTLY SECURED TO A STANDARDISED FALL ARREST SYSTEM.

In case of device failure or malfunction, it must be replaced immediately with equivalent equipment, thus ensuring continuity of emergency procedures.

The device must be maintained in accordance with the manufacturer's instructions, including periodic servicing, recommended maintenance and adherence to the estimated service life specified in the User's Manual.

If you have any doubts about the correct interpretation of the manual or about the use of the device, it is recommended that you contact IRUDEK to obtain the necessary clarifications.

It is essential to ensure that all elements of the system are operating properly before using the equipment. Under no circumstances should a configuration be used in which the operation of one component interferes with or is affected by another.

All components of the device must comply with the applicable regulations and the specific instructions for use. The device may only be used for the purposes intended by the manufacturer.

The use of the equipment is prohibited for persons whose health condition may compromise their safety during daily use or in an emergency situation.

It is mandatory to have an emergency action plan, properly prepared and available for immediate implementation in case of need.

Any technical intervention on the equipment, whether repair or modification, may only be carried out by the manufacturer or an authorised representative. It is expressly forbidden to make any adjustments or alterations without the written consent of the manufacturer.

THE EMAN AND EMAN ROLL STRETCHER IS NOT DESIGNED AS A SPINAL IMMOBILISATION DEVICE. IF THE CASUALTY HAS TRAUMA CONSISTENT WITH SPINAL INJURIES, IMMOBILISATION WITH A SPINAL IMMOBILISATION DEVICE IS APPROVED FOR THIS FUNCTION MUST BE ENSURED BEFORE THE CASUALTY IS PLACED ON THE EMAN AND EMAN ROLL STRETCHER.

Rescue procedures using ropes, harnesses, spider harnesses and other specialised equipment require specific skills and involve risks. Rescue personnel should not perform them without prior professional guidance and training.

Ensure that the person to be rescued is correctly positioned so that all rescue manoeuvres can be carried out safely. At all times, care must be taken to prevent the webbing straps from slipping and creating a dangerous situation in the event of an unexpected movement (e.g. unhooking at close range).

For safety reasons, and before using the equipment, check that there are no obstacles that could interfere with the normal course of the intervention. Since this type of manoeuvre is usually carried out in emergency situations, it is essential that rescuers are properly trained and have previously practised this type of operation.

The safety of the user depends on the correct operation of the equipment and proper understanding of the instructions in this manual.

This device is designed for outdoor use on construction sites. It may only be operated by trained and fit persons under the supervision of qualified personnel.

Certain medical conditions may affect the safety of the user; if in doubt, consult a physician.

Avoid using the equipment in contact with sharp edges, small diameter structures or areas susceptible to corrosion, as this may compromise the strength of the belts.

The addition or replacement of any component is prohibited.

When using the rescue loop EN 1498 included with the EMAN ROLL stretcher independently of the stretcher Fig. 23, remember to consider the risk of trauma caused by being suspended in a vacuum; therefore, the time the person is restrained with the rescue loop should be as short as possible, in order to reduce the effects on breathing and circulation.

DISINFECTION AND CLEANING OF THE STRETCHER

To reduce the risk of contracting blood-borne diseases, follow the cleaning and disinfection instructions in this manual when using the stretcher.

Rescue procedures using ropes, harnesses, spider harnesses and other specialised equipment require specific skills and involve risks. Rescue personnel should not perform them without prior professional guidance and training.

For plastic and metal parts: Apply a disinfectant cleaner to all parts of the stretcher, following the instructions on the cleaner. Before storing, make sure all parts are completely dry.

For textile parts: Close the belts before cleaning to prevent damage in washing machines. Place the belts in a mesh bag (not included) and place them in an industrial washing machine. Wash them in warm water, using a disinfectant soap and selecting the delicate cycle. Remove the belts from the mesh bag and allow them to air dry. Do not use tumble dryers.

LIFELINE EN 1498

The rescue lanyard EN 1498:2006 class A is a lanyard with sewn terminals that should only be used to rescue conscious adults and evacuate them upwards or downwards with the intervention of a professional who accompanies the casualty at all times.

It can be used independently or as a strap to be placed on the upper part of the stretcher, as a suspension strap placed at shoulder height.

NOMENCLATURE

10. Rescue loop, 11. lanyard, 12. anchorage points

USE

Pass the lifeline loop under the person's arm and tighten it around the torso. Check that the strap is not twisted and that it is correctly positioned on the back. Finally, connect the end of the lanyard to the evacuation system using an EN392 compliant connector. Fig. 23

The maximum permitted load is 150kg.

The anchorage point of rescue systems must offer sufficient strength; the minimum required strength is 12 kN (EN 952:2012). Ensure that manoeuvres are carried out in such a way as to limit the pendulum effect as well as the risk of falling.

For use on the stretcher, please refer to the instructions for use of the stretcher.

COMPATIBILITIES

This rescue loop is designed to be used only with compatible evacuation systems (e.g. those conforming to EN341 and EN1496). Creating a customised system can be dangerous, as one safety function may interfere with another. For this reason, it is important to consult the instructions for use for each component of the system before use.

MATERIAL

- Material of the tapes: polyester

The rescuer must ensure that the person being rescued is not exposed to hazards caused by displacement of the rescue loop straps or contact with connecting elements.

The rescued person should remain attached to the device for as short a time as possible to reduce the risk associated with breathing and circulation. Remaining suspended for a prolonged period of time can result in serious impaired circulation and breathing. It is essential to act quickly and follow proper rescue procedures.

CHECKS BEFORE USE

Prior to use, a visual and functional inspection of its components must be carried out by the user, verifying that they do not show signs of deterioration, excessive wear, corrosion, abrasions, degradation due to UV radiation, cuts and incorrect use. Special attention should be paid to straps, seams, anchorage rings, buckles and adjustment elements.

If any defects, anomalies or damage are found in the personal protection equipment that entail a loss of protection, it must be removed from use.

CHECKS DURING USE

While using the equipment, pay special attention to any hazardous circumstances that may affect equipment performance and user safety, including the following:

- The labelling on the safety components.
- Accidental contact with sharp edges.
- Various types of damage, such as cuts, abrasion and/or corrosion.
- The negative effect of weather conditions.
- "Pendulum" falls.
- Effects of extreme temperatures.
- Effects after contact with chemical products.
- Electrical conductivity.
- It is essential that all fasteners and fittings are checked regularly.

WARRANTY

This product has a 3-year warranty that covers manufacturing and raw material defects. The warranty does not cover wear, corrosion or damage caused by storage, transport or improper or intensive use.

The warranty application must be submitted along with the purchase receipt. If a manufacturing defect is found, IRUDEK agrees to repair, replace or refund the product for an amount that does not exceed the price stated in the product invoice.

WASTE MANAGEMENT

Products without electrical components: dispose of the product safely at the end of its useful life. Separate textiles, plastics and metal materials as far as possible for environmental management.

Electrical or electronic products / with batteries: This product contains electrical components or batteries and must not be disposed of with household waste. Please hand it over to an authorised waste collector or consult www.irudek.com for proper disposal.

USEFUL LIFE

The estimated useful life of textile equipment is 12 years from the date of manufacture (2 years of storage and 10 years of use). Metal equipment has an unlimited useful life.

The following factors can reduce the product's useful life: intensive use, contact with chemical substances, especially aggressive environments, exposure to extreme temperatures, exposure to ultraviolet rays, abrasion, cuts, strong impacts, improper use, transport and/or maintenance.

TRANSPORT

This personal protection equipment must be transported in packaging that protects it against humidity and any mechanical, chemical and/or thermal damage.

STORAGE

This personal protection system must be stored in a package with plenty of room in a dry place, protected against sunlight, ultraviolet rays, dust, sharp objects, extreme temperatures and aggressive substances.

REQUIREMENTS

Prior to the use of the equipment, a rescue plan has to be established in order to be able to execute it in case of emergency.

Do not make any changes or add any elements to the equipment without prior written authorisation from the manufacturer.

The equipment must not be used outside its scope of limitations or for any purpose other than its intended purpose.

Make sure that the equipment components are compatible with the system it is assembled to. Make sure that the elements are appropriate for the proposed application. It is forbidden to use the protection system if the operation of an individual component is affected by or interferes with the operation of another component. Perform a periodic inspection of the connections and adjustments of the components to ensure that they do not come loose accidentally.

If any wear or damage is detected or there are any doubts as to safe conditions of use, this personal protection equipment should be removed from use immediately. It must not be used again until an authorised individual presents a written confirmation that it is in suitable condition to be used.

If the equipment has prevented a fall, it should be removed from service.

Before each use, for safety purposes it is essential to verify the minimum distance of free space required under the user's feet to avoid colliding with the ground or any other obstacle in the event of a fall. Detailed information regarding the minimum requirements of free space can be found in the instructions of the corresponding fall prevention system components.

If the product is resold outside the original country of destination, the reseller must provide instructions of use, maintenance, periodic inspection and repair in the language of the country where the equipment will be used.

MAINTENANCE INSTRUCTIONS

Visual inspection

Users should perform a visual and functional inspection of the equipment before using it.

If the equipment has undergone unusual or extraordinary conditions, a special inspection should be carried out by the manufacturer or a competent person authorised by the manufacturer.

At least every 12 months, a thorough periodic overhaul must be carried out by the manufacturer or a competent person authorised by the manufacturer, in strict accordance with IRUDEK's periodic overhaul procedures. The safety of the users depends on the continued efficiency and durability of the equipment. The periodic inspection must be certified according to the requirements of EN365:2004, determining the validity of the certificate and the date of the next inspection.

The product marking must be legible.

Any pertinent observations must be entered in the equipment inspection certificate.

If any defects, anomalies or damage are found in the personal protection equipment that entail a loss of protection, it must be removed from use.

Cleaning

This personal protection equipment must be cleaned without causing any damage to the materials used for its manufacture or to its use. The cleaning procedure must be followed strictly. Clean textile and plastic materials (belts, ropes) with a cotton or cloth or a brush. Do not use any type of abrasive material. To clean the equipment thoroughly, wash it by hand at a temperature between 30 and 40°C, using neutral soap. Use a moist cloth for the metal parts. If the equipment gets wet due to use or cleaning, let it dry naturally in a well-ventilated place, away from direct heat or chemical compounds.

The disinfection process shall be carried out in the same way as the deep cleaning process.

Repair

The equipment must only be repaired by the manufacturer or a person authorised to do so and following the procedures established by the manufacturer. Instructions for repair will be provided in the official languages of the country where the equipment is put to use.

CONTROL SHEET

The control sheet should be completed before the equipment is delivered for its first use.

All the information about the personal protection equipment (name, serial number, date of purchase and date of first use, user name, periodic inspection and repair log and next periodic inspection date) must be entered in the equipment's control sheet.

The sheet must be completed exclusively by the person responsible for the protection equipment.

Iruckek

The Iruckek application is used for easy, effective control of fall prevention equipment. Its use is recommended to trace these products, thereby replacing the Control Sheet.

TRANSLATIONS: EXPLANATORY NOTE

The translation of all documents originally written in Spanish is done by an external translator and is provided as part of an information service to the global community. Inaccuracies may arise as a result of language restrictions and translation errors. IRUDEK does not verify the accuracy of translations made by third parties and therefore assumes no liability whatsoever in relation to any disputes and/or claims that may arise as a result of errors, omissions or ambiguities in the translated material contained herein. Any person or body relying on such translated material does so at his or her own risk and responsibility. In case of doubt or dispute as to the accuracy of the translated text, the original language equivalent shall prevail. If you wish to report an error or inaccuracy in the translation, we invite you to write to us at info@irudek.com

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