
IRUDEK

Manual de uso ES
User manual EN
Manuale d'uso IT
Manual do utilizador PT
Manual do utilizador BR
Benutzerhandbuch DE
Manuel de l'utilisateur FR
Podręcznik użytkownika PL
Felhasználói kézikönyv HU
Manual de utilizare RO
Användarmanual SV
Brukerhåndbok NO
Käyttöohje FI
Naudotojo vadovas LT
Kasutusjuhend ET
Ръководство за потребителя BG
Упутство за употребу SR
Používateľská príručka SK
Kullanım Kılavuzu TR

IRUDEK 2000 S.L.
Pol. Erribera 8A
20150 Aduna (Guipúzcoa)
España
Tfno: +34 943 69 26 17
Fax: +34 943 69 25 26
irudek@irudek.com

NOMENCLATURE:


Fig.1


Fig.2


Fig.3


Fig.4


Fig.5

Fig.6


Fig.7


Fig.8


Fig.9


Fig.10


Fig.11


Fig.12


Fig.13

Fig.14

Fig.15


Fig.16

Fig.17

Fig.18


Fig.19


Fig.22

Fig.20


Fig.23




Fig.21

MAINTENANCE

CONTROL SHEET

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EN

Read the operating instructions carefully before using the stretcher, train yourself properly, familiarise yourself with it and use it responsibly. Activities at height involve serious risks not outlined in this manual, where each user is responsible for the management of such risks, their safety, their actions and the consequences of these, if you do not assume this or do not understand this manual, do not use the equipment.

MEDICAL CONDITIONS THAT MAY AFFECT THE SAFETY OF THE USER:

UNDER NORMAL CONDITIONS OF USE, CARDIOVASCULAR PROBLEMS, RESPIRATORY DISEASES, MUSCULOSKELETAL DISORDERS AFFECTING THE SPINE, HIPS OR KNEES, OBESITY OR EXCESS WEIGHT, AS WELL AS NEUROLOGICAL OR BALANCE DISORDERS AND PERIPHERAL CIRCULATORY PROBLEMS HINDERING VENOUS RETURN MUST BE TAKEN INTO ACCOUNT. IN AN EMERGENCY SITUATION FOLLOWING THE ARREST OF A FALL, THE RISKS ARE PARTICULARLY SERIOUS IN PEOPLE PREDISPOSED TO HARNESS SYNDROME, WITH A HISTORY OF TRAUMA OR INJURY, WITH COAGULATION PROBLEMS OR ON ANTICOAGULANT TREATMENT, AND IN THOSE WITH NEUROLOGICAL CONDITIONS THAT CAN CAUSE SEIZURES.

WARNINGS

Improper use of the Stretcher may result in accidents or injury. Use it only for the purposes indicated in this manual and according to the instructions provided.

Only properly trained personnel should handle the stretcher. Use by untrained persons may result in injury to themselves and the patient.

During the transfer, the patient must be accompanied at all times, as leaving the patient unsupervised can be dangerous.

Always secure the patient with the appropriate belts. An unrestrained patient could fall from the stretcher and suffer injury. Use a restraint system appropriate to the mission and keep the operation of the stretcher under control, directing the assistants correctly to avoid accidents.

Improper maintenance of the Stretcher can also lead to hazards. Only carry out maintenance in accordance with the instructions in this manual.

The use of unapproved parts, services or accessories can be dangerous. Use only IRUDEK approved components and services, as well as accessories certified by IRUDEK for use with the Stretcher.

DESCRIPTION

The EMAN stretcher is a nest stretcher that allows the casualty to be transported both horizontally and vertically.

For vertical manoeuvres, the casualty must have a fall arrest system independent of the stretcher.

The stretchers are made of highly durable material and have dimensions of 200 x 90 cm.

It is possible to separate the stretcher into two pieces, in a way that facilitates its transport.

Height (mm)	200
Width (mm)	650
Length (mm)	2140
Split length (mm)	1180
Weight (without accessories) (kg)	kg 13.8
Load capacity (kg)	200

Content of the set:

4 carabiners, 4 straps for closing the stretcher, 1 suspension system, 1 stirrup and its anchoring strap, 1 perimeter cord, 1 upper body, 1 lower body, 2 side locking pins, 1 lower locking pin, 1 lower locking pin, 1 upper body, 1 lower body, 1 upper body, 1 lower body, 2 side locking pins, 1 lower locking pin, 1 lower locking pin.

Information on this equipment is available at the following link: http://www.irudek.com

NOMENCLATURE

1. Suspension straps 2. Stretcher body, 3. Perimeter rope, 4. Eyelets for suspension straps, 5. Handles, 6. Footrest, 7. Fastening straps, 8. Comfort mattress, 9. Connecting pins, 10. Connecting eyelets, 11. Locking pins

ASSEMBLY

Fig.1: Remove the stretcher from the box.

Fig.2: Insert the white plastic tubes into the corresponding receiving hole.

Fig.3: Check that the assembly between the two parts has been carried out correctly.

Fig.4 After inserting the stainless steel pins shown in the circle into the corresponding holes.

Fig.5: On the back side you can see the pins going through the hole.

Fig.6: Then insert the buckle attached to the pin into the small hole in the pin to complete the attachment.

Fig.7: On the underside of the basket stretcher there is also a fixing device.

Fig. 8: Insert the large pin, as shown in the figure, into the hole marked by the circle.

Fig.9-10: Insert the latch into the pin hole to complete the fixing.

Fig.11: To install the footrest, thread the orange tape through the eyelets of the plastic part.

Fig.12: Insert the metal anchor into the side hole of the stretcher.

Fig.13: The footrest attachment is adjustable and should be as shown in the picture.

Fig.14-15: For the installation of the casualty restraint straps, pass the free end of the strap through the perimeter rope.

Fig.16: Place the free end of the strap through the locking buckle as shown in the picture.

Fig.17: Adjust the length of the tape.

Fig.18: Once in place, they should be in the position shown.

Fig.19: To install the lifting straps, attach the carabiner to the thin strap attached to the buckle.

Fig.20 the system should look as shown in the picture.

Fig.21: Pass the end of the wide strap through the metal eyelets of the stretcher.

Fig.22: Attach the end of the wide strap to the carabiner.

Fig.23: Check that the straps are not crossed or twisted.

RESTRICTIONS ON USE

The device must only be operated by suitably trained personnel.

THE STRETCHER IS NOT A FALL-ARREST SYSTEM; THEREFORE, WHEN CARRYING OUT MANOEUVRES INVOLVING A RISK OF FALLING BETWEEN DIFFERENT LEVELS, THE PERSON ON THE STRETCHER MUST BE CORRECTLY SECURED TO A STANDARD FALL-ARREST SYSTEM.

In the event of a fault or malfunction, the equipment must be replaced immediately with an equivalent unit, thereby ensuring the continuity of emergency procedures.

The device must be maintained in accordance with the manufacturer's instructions, including regular inspections, recommended maintenance and adherence to the estimated service life specified in the User Manual.

If you have any queries regarding the correct interpretation of the manual or the use of the device, we recommend that you contact IRUDEK for the necessary clarification.

It is essential to ensure that all system components are functioning correctly before using the equipment. Under no circumstances should a configuration be used in which the operation of one component interferes with or is affected by another.

All components of the equipment must comply with current regulations, as well as with the specific instructions for use. The equipment must only be used for the purposes intended by the manufacturer.

Use of the equipment is prohibited for people whose health condition may compromise their safety during everyday use or in an emergency.

It is compulsory to have an emergency action plan in place, which must be properly prepared and available for immediate implementation should the need arise.

Any technical work on the equipment, whether repairs or modifications, may only be carried out by the manufacturer or an authorised representative. It is expressly forbidden to make any adjustments or alterations without the manufacturer's written consent.

THE EMAN AND EMAN ROLL STRETCHERS ARE NOT DESIGNED AS SPINAL IMMOBILISATION DEVICES. IF THE CASUALTY HAS INJURIES CONSISTENT WITH SPINAL TRAUMA, THEY MUST BE IMMOBILISED USING A SPINAL IMMOBILISATION DEVICE APPROVED FOR THAT PURPOSE BEFORE BEING PLACED ON THE EMAN AND EMAN ROLL STRETCHERS.

Rescue procedures involving ropes, harnesses, spider harnesses and other specialist equipment require specific skills and involve risks. Rescue personnel must not carry out such procedures without first receiving professional guidance and training.

Check that the person to be rescued is correctly positioned so that all rescue manoeuvres can be carried out safely. Care must be taken at all times to ensure that the straps do not slip, as this could create a hazardous situation in the event of an unexpected movement (for example, the harness coming undone at close range).

For safety reasons, and before using the equipment, check that there are no obstacles that could interfere with the normal course of the operation. As this type of manoeuvre is usually carried out in emergency situations, it is essential that rescuers are properly trained and have practised this type of operation beforehand.

User safety depends on the equipment functioning correctly and on a proper understanding of the instructions in this manual.

This device is designed for outdoor use on construction sites. It must only be operated by trained personnel who are physically fit and under the supervision of equally qualified staff.

Certain medical conditions may affect the user's safety; if in doubt, consult a doctor.

Avoid using the equipment in contact with sharp edges, structures with a small diameter or areas susceptible to corrosion, as these may compromise the strength of the belts.

It is prohibited to add or replace any component.

DISINFECTION AND CLEANING OF THE TREATMENT TABLE

To reduce the risk of contracting diseases transmitted through contact with blood, follow the cleaning and disinfection instructions in this manual when using the treatment table.

Rescue procedures involving ropes, harnesses, spider harnesses and other specialist equipment require specific skills and involve risks. Rescue personnel must not carry out such procedures without first receiving professional guidance and training.

For plastic and metal parts: Apply a disinfectant cleaner to all parts of the treatment table, following the cleaner's instructions. Before storing it, ensure that all parts are completely dry.

For textile parts: Fasten the belts before cleaning them to prevent them from being damaged in the washing machines. Place the belts in a mesh bag (not included) and put them in an industrial washing machine. Wash them in lukewarm water, using a disinfectant soap and selecting the delicate cycle. Remove the belts from the mesh bag and leave them to air dry. Do not use tumble dryers.

CHECKS BEFORE USE

Prior to use, a visual and functional inspection of its components must be carried out by the user, verifying that they do not show signs of deterioration, excessive wear, corrosion, abrasions, degradation due to UV radiation, cuts and incorrect use. Special attention should be paid to straps, seams, anchorage rings, buckles and adjustment elements.

If any defects, anomalies or damage are found in the personal protection equipment that entail a loss of protection, it must be removed from use.

CHECKS DURING USE

While using the equipment, pay special attention to any hazardous circumstances that may affect equipment performance and user safety, including the following:

- The labelling on the safety components.
- Accidental contact with sharp edges.
- Various types of damage, such as cuts, abrasion and/or corrosion.
- The negative effect of weather conditions.
- "Pendulum" falls.
- Effects of extreme temperatures.
- Effects after contact with chemical products.
- Electrical conductivity.
- It is essential that all fasteners and fittings are checked regularly.

WARRANTY

This product has a 3-year warranty that covers manufacturing and raw material defects. The warranty does not cover wear, corrosion or damage caused by storage, transport or improper or intensive use.

The warranty application must be submitted along with the purchase receipt. If a manufacturing defect is found, IRUDEK agrees to repair, replace or refund the product for an amount that does not exceed the price stated in the product invoice.

WASTE MANAGEMENT

Products without electrical components: dispose of the product safely at the end of its useful life. Separate textiles, plastics and metal materials as far as possible for environmental management.

Electrical or electronic products / with batteries: This product contains electrical components or batteries and must not be disposed of with household waste. Please hand it over to an authorised waste collector or consult www.irudek.com for proper disposal.

USEFUL LIFE

The estimated useful life of textile equipment is 12 years from the date of manufacture (2 years of storage and 10 years of use). Metal equipment has an unlimited useful life.

The following factors can reduce the product's useful life: intensive use, contact with chemical substances, especially aggressive environments, exposure to extreme temperatures, exposure to ultraviolet rays, abrasion, cuts, strong impacts, improper use, transport and/or maintenance.

TRANSPORT

This personal protection equipment must be transported in packaging that protects it against humidity and any mechanical, chemical and/or thermal damage.

STORAGE

This personal protection system must be stored in a package with plenty of room in a dry place, protected against sunlight, ultraviolet rays, dust, sharp objects, extreme temperatures and aggressive substances.

REQUIREMENTS

Prior to the use of the equipment, a rescue plan has to be established in order to be able to execute it in case of emergency.

Do not make any changes or add any elements to the equipment without prior written authorisation from the manufacturer.

The equipment must not be used outside its scope of limitations or for any purpose other than its intended purpose.

Make sure that the equipment components are compatible with the system it is assembled to. Make sure that all the elements are appropriate for the proposed application. It is forbidden to use the protection system if the operation of an individual component is affected by or interferes with the operation of another component. Perform a periodic inspection of the connections and adjustments of the components to ensure that they do not come loose accidentally.

If any wear or damage is detected or there are any doubts as to safe conditions of use, this personal protection equipment should be removed from use immediately. It must not be used again until an authorised individual presents a written confirmation that it is in suitable condition to be used.

If the equipment has prevented a fall, it should be removed from service.

Before each use, for safety purposes it is essential to verify the minimum distance of free space required under the user's feet to avoid colliding with the ground or any other obstacle in the event of a fall. Detailed information regarding the minimum requirements of free space can be found in the instructions of the corresponding fall prevention system components.

If the product is resold outside the original country of destination, the reseller must provide instructions of use, maintenance, periodic inspection and repair in the language of the country where the equipment will be used.

MAINTENANCE INSTRUCTIONS

Visual inspection

Users should perform a visual and functional inspection of the equipment before using it.

If the equipment has undergone unusual or extraordinary conditions, a special inspection should be carried out by the manufacturer or a competent person authorised by the manufacturer.

At least every 12 months, a thorough periodic overhaul must be carried out by the manufacturer or a competent person authorised by the manufacturer, in strict accordance with IRUDEK's periodic overhaul procedures. The safety of the users depends on the continued efficiency and durability of the equipment. The periodic inspection must be certified according to the requirements of EN365:2004, determining the validity of the certificate and the date of the next inspection.

The product marking must be legible.

Any pertinent observations must be entered in the equipment inspection certificate.

If any defects, anomalies or damage are found in the personal protection equipment that entail a loss of protection, it must be removed from use.

Cleaning

This personal protection equipment must be cleaned without causing any damage to the materials used for its manufacture or to the user. The cleaning procedure must be followed strictly. Clean textile and plastic material's (belts, ropes) with a cotton or cloth or a brush. Do not use any type of abrasive material. To clean the equipment thoroughly, wash it by hand at a temperature between 30 and 40ºC, using neutral soap. Use a moist cloth for the metal parts. If the equipment gets wet due to use or cleaning, let it dry naturally in a well-ventilated place, away from direct heat or chemical compounds.

The disinfection process shall be carried out in the same way as the deep cleaning process.

Repair

The equipment must only be repaired by the manufacturer or a person authorised to do so and following the procedures established by the manufacturer. Instructions for repair will be provided in the official languages of the country where the equipment is put to use.

CONTROL SHEET

The control sheet should be completed before the equipment is delivered for its first use.

All the information about the personal protection equipment (name, serial number, date of purchase and date of first use, user name, periodic inspection and repair log and next periodic inspection date) must be entered in the equipment's control sheet.

The sheet must be completed exclusively by the person responsible for the protection equipment.

IruCheck

The IruCheck application is used for easy, effective control of fall prevention equipment. Its use is recommended to trace these products, thereby replacing the Control Sheet.

TRANSLATIONS: EXPLANATORY NOTE

The translation of all documents originally written in Spanish is done by an external translator and is provided as part of an information service to the global community. Inaccuracies may arise as a result of language restrictions and translation errors. IRUDEK does not verify the accuracy of translations made by third parties and therefore assumes no liability whatsoever in relation to any disputes and/or claims that may arise as a result of errors, omissions or ambiguities in the translated material contained herein. Any person or body relying on such translated material does so at his or her own risk and responsibility. In case of doubt or dispute as to the accuracy of the translated text, the English language equivalent shall prevail. If you wish to report an error or inaccuracy in the translation, we invite you to write to us at info@irudek.com

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